



Protexer

Personal Protection Apparel
Protection. Comfort. Infection Prevention.

TT TEAM
Medical
A DIVISION OF **TT TEAM** Technologies

PROTE~~X~~ER

Protexer | PPE Disposable Gown

SMS – Isolation Gown: ICP-8150, ICP-8300, ICP-8150FC, RR-200, RR-203XL

Recommended for:
House wide/Med-Surg/ICU

- Easy to don over the head gown
- Thumb loop with elastic cuff offers double protection at the wrist to ensure no gap between the glove and gown.
- Universal & 3XL Sizes
- Choose between Open Back or Full Coverage
- Side Ties on full coverage gown
- Constructed from medium weight material
- Meets AAMI PB70:2012 Level-II standards for full coverage



Procedure Dental Gown: ICP-8200, ICP-8250

Recommended for:
Dental Clinic and Dental Procedure

- Full Coverage wrap gown
- Universal and 2XL Sizes
- Comfortable Knit Cuff
- Velcro neck closure
- Constructed from medium weight SMS material
- Meets AAMI PB70:2012 Level-II standards



Protexer | PPE Disposable Gown

Poly-coated Isolation Gown: ICP-5050Y, ICP-5050YFC, RR-300

Recommended for:
Isolation/ICU/Emergency/Maternity

- Easy to don over the head gown
- Thumb loop with elastic cuff offers double protection at the wrist to ensure no gap between the glove and gown.
- Universal Size
- Choose between Open Back or Full Coverage
- Constructed from Polycoated material (PPSB+PE coating)
- Polyethylene coated materials tested against ASTM F1670



Open Back



Full Coverage

SMS – Isolation Gown: RR-300FC

Recommended for:
House wide/Med-Surg/ICU

- Easy to don over the head gown
- Thumb loop with elastic cuff offers double protection at the wrist to ensure no gap between the glove and gown.
- Universal
- Full Coverage
- Constructed from heavy weight material
- Meets AAMI PB70:2012 Level-III standards for full coverage



Full Coverage

Protexer | PPE Disposable Gown

Procedure Gown Tested for Chemotherapy Drugs: ICP-7050C

Recommended for:
Chemotherapy/ICU/Endoscopy/GI Lab

- Easy to don over the head gown
- Universal Size
- Comfortable Knit Cuff
- Open Back
- Constructed from polycoated material (PPSB+PE coating)
- Polyethylene coated materials tested against ASTM F1670
- Chemotherapy drug tested according to ASTM F739-12 with no breakthrough observed during 360 minutes of testing
 - Doxorubicin Hydrochloride
 - Drug Tested
 - Etoposide (Toposar)
 - Cyclophosphamide (Cytosan)
 - Fluorouracil
 - Paclitaxel (Taxol)
 - Thiotepa
 - Vincristine Sulfate
 - Carmustine (BCNU)
 - Cisplatin
 - Dacarbazine (DTIC)



Open Back

Protexer | PPE Disposable Gown

Procedure Gown: ICP-7052FCC, ICP-7053FCC

Recommended for:
ICU/Endoscopy/GI Lab

- Easy to don over the head gown
- XL & Universal Sizes
- Comfortable Knit Cuff
- Coverage: Full Coverage
- Constructed from polycoated material (PPSB+PE coating)
- Polyethylene coated material tested against ASTM F1670



Polyethylene – Isolation Gown: ICP-3845

Recommended for:
**Decontamination/Emergency/
ICU/Isolation**

- Easy to don over the head gown
- Universal Size
- Thumb loop at the wrist to ensure no gap between the glove and gown
- Coverage: Open Back
- Constructed from Polyethylene material



Protexer | PPE Disposable Gown

Protexer Gown Specifications

VA/DoD Item #	Description	AAMI Level	Material	Color	Size	Coverage	Cuff	Neck	Waist	Latex Free	Packaging String	Trade Approved Country of Origin
ICP-8150FC	SMS Isolation Gown, Yellow, Universal, Full Coverage, AAMI 2	II	SMS	Yellow	Universal	Full Coverage	Thumbloop w/elastic wrist	Over the head	Side Tie	Yes	10ea/10pk/100cs	X
ICP-8150	SMS Isolation Gown, Yellow, Universal Size, Open Back		SMS	Yellow	Universal	Open Back	Thumbloop w/elastic wrist	Over the head	Apron Tie	Yes	10ea/10pk/100cs	X
RR-200	SMS Isolation Gown, Yellow, Universal Size, Open Back		SMS	Yellow	Universal	Open Back	Thumbloop w/elastic wrist	Over the head	Apron Tie	Yes	10ea/10pk/100cs	
ICP-8300	SMS Isolation Gown, Yellow, 3XL, Open Back		SMS	Yellow	3XL	Open Back	Thumbloop w/elastic wrist	Over the head	Apron Tie	Yes	10ea/10pk/100cs	X
RR-203XL	SMS Isolation Gown, Yellow, 3XL, Open Back		SMS	Yellow	3XL	Open Back	Thumbloop w/elastic wrist	Over the head	Apron Tie	Yes	10ea/10pk/100cs	
ICP-5050Y	Polycoated Isolation Gown, Yellow, Universal Size, Open Back		Polycoated	Yellow	Universal	Open Back	Thumbloop w/elastic wrist	Over the head	Apron Tie	Yes	10ea/10pk/100cs	X
RR-300	Polycoated Isolation Gown, Yellow, Universal Size, Open Back		Polycoated	Yellow	Universal	Open Back	Thumbloop w/elastic wrist	Over the head	Apron Tie	Yes	10ea/10pk/100cs	
ICP-5050YFC	Polycoated Isolation Gown, Yellow, Universal, Full Coverage	II	Polycoated	Yellow	Universal	Full Coverage	Thumbloop w/elastic wrist	Over the head	Side Tie	Yes	10ea/10pk/100cs	X
RR-300FC	SMS Isolation Gown, Yellow, Universal, Full Coverage	III	SMS	Yellow	Universal	Full Coverage	Thumbloop w/elastic wrist	Over the head	Side Tie	Yes	10ea/10pk/100cs	
ICP-3845	Polyethylene Gown, Blue, Universal Size, Open Back		Polyethylene	Blue	Universal	Open Back	Thumbloop w/elastic wrist	Over the head	Apron Tie	Yes	15ea/5bx/75cs	X
ICP-7050C	Procedure Gown, White, Knit Cuffs, Universal Size, Open Back, Tested for Chemo Drugs		PP/PE	White	Universal	Open Back	Knit Cuffs	Over the head	Apron Tie	Yes	25ea/2pk/50cs	X
ICP-7052FCC	Procedure Gown, White, Knit Cuffs, Extra Large, Full Coverage		PP/PE	White	XL	Full Coverage	Knit Cuffs	Over the head	Back Tie	Yes	25ea/2pk/50cs	X

*PB70:2012: AATCC 42-Spray Impact Penetration Test and AATCC 127 Hydrostatic Pressure Test

*Latex Free - Not made with natural rubber

Protexer | PPE Disposable Gown

Protexer Gown Specifications (cont'd)

VA/DoD Item #	Description	AAMI Level	Material	Color	Size	Coverage	Cuff	Neck	Waist	Latex Free	Packaging String	Trade Approved Country of Origin
ICP-7053FCC	Procedure Gown, White, Knit Cuffs, Universal, Full Coverage		PP/PE	White	Universal	Full Coverage	Knit Cuffs	Over the head	Back Tie	Yes	25ea/2pk/50cs	X
ICP-8200	Dental Gown, Yellow, Full Coverage, Velcro Neck, Knit Cuffs, Universal	II	SMS	Yellow	Universal	Full Coverage	Knit Cuffs	Velcro Neck	Side Tie	Yes	10ea/10pk/100cs	X
ICP-8250	Dental Gown, Yellow, Full Coverage, Velcro Neck, Knit Cuffs, 3XL	II	SMS	Yellow	2XL	Full Coverage	Knit Cuffs	Velcro Neck	Side Tie	Yes	10ea/10pk/100cs	X

*PB70:2012: AATCC 42-Spray Impact Penetration Test and AATCC 127 Hydrostatic Pressure Test

*Latex Free - Not made with natural rubber

Protexer Gowns Offer Superior Comfort, Rapid Release and Other Features for Increased Staff Compliance

Simplicity: Universal size fits most employees. Unique folding and packaging maximize supply and minimizes restocking expense.

Protection: Elastic wrists and thumb loops ensure complete, secure coverage.

Style: Available in SMS or Poly-coated materials with open back or full coverage design

Convenience: Over-the-head donning makes compliance easy and perforations at the shoulders simplify removal and minimize aerosolization.

Protexer | PPE Disposable Gown

Effectively Donning and Removing Protexer Gowns

Helpful Hints for Donning

1. Locate label at the top of the gown. Hold at neck and allow to unfold.
The patented fold allows the gown to unfold in the ideal position for donning.
2. Place head through apron style opening.
3. With thumbs facing upward, slide one arm at a time into sleeves through the thumb loop design.
4. Pull the right side of the gown behind your back and then overlap the left.
5. Gently secure ties behind the back.

Please note that any article worn on the arms, wrists or hands such as a watch or ring, could tear the gown. Please remove all adornments before donning the gown.

Helpful Hints for Removing

1. To remove, locate the base of the tie straps at hip and pull to release one at a time.
(Do NOT pull both sides at once as you risk aerosolization of contaminants on the gown).
2. Reach up to the shoulder, grasp the gown and pull down.
Repeat on the other side. (Do NOT pull both sides at once as you risk aerosolization of contaminants on the gown).
3. Roll gown into itself removing gloves into the gown being careful to only touch the inside of the gown once the gloves are removed.
4. Properly dispose of the gown.



Protexer | PPE Disposable Comfort Scrubs

Disposable Protexer Scrubs

Disposable Protexer Scrubs are a fresh alternative to laundered scrubs. The breathable fabric offers supreme comfort and the triple layer fluid-resistance maximizes protection. Whether as an alternative to laundered scrubs or as an emergency back-up supply, consider Protexer Scrubs for those inevitable occasions when theft/shrinkage, excess wear, or a surge depletes conventional inventory.

- Top: V-neck and Breast Pocket
- Pants: Elastic Waist
- Fluid resistant
- 3-layer SMS, breathable fabric

Adult Protexer Scrubs

ICP-9000TOP Comfort Scrub Top, Regular, Blue, 100/Case

ICP-9000BTM Comfort Scrub Bottom, Regular, Blue, 50/Case

ICP-9050TOP Comfort Scrub Top, Extra Large, Blue, 100/Case

ICP-9050BTM Comfort Scrub Bottom, Extra Large, Blue, 50/Case

Pediatric Scrubs

- Top: V-neck
- Pants: Drawstring is sewn in place in the back. Breaks with 20 pounds of pressure.

RR-4000S Small, 50 sets/case

RR-4000M Medium, 50 sets/case

RR-4000L Large, 50 sets/case

RR-4000XL X-large, 50 sets/case



Comfort Scrubs



Pediatric Scrubs

Protexer | Lab Coats

Lab Coats

Protexer Lab Coats maximize safety and compliance for both patients and healthcare workers. Our Lab Coats help reduce potential for cross-contamination while offering superior protection and excellent comfort.

- 3-layer SMS fluid resistant fabric
- 3 Pockets
- 6 Snaps
- Color: Blue
- Knit Cuff
- Traditional Collar
- 50 per case

ICP-9054LC Universal

ICP-9055LC 2XL

ICP-9057LC 4XL



Protexer | Disposable Personal Protection Equipment

Disposable Shoe Covers

ICP-SCNS0001, ICP-SCNS0002, ICP-SCNS0003

The ICP disposable shoe covers are made of a polypropylene fabric and are 100% latex and silicone free. They provide great fluid protection and slip resistances.

ICP-SCNS0001 & ICP-SCNS0003 Features

- Elastic band for secure fit
- Lightweight but durable
- Non-Skid pattern
- Large and Universal Sizes
- Fluid resistant
- Color: Blue
- 300/case - ICP-SCNS0001
- 200/case - ICP-SCNS0003

ICP-SCXT0002 Features

- Elastic band for secure fit
- Lightweight but durable
- Extra traction Non-Skid pattern
- XL Sizes
- Fluid resistant
- Color: Blue
- 200/case



PR-SCNS03 & PR-SCNS04

The PR-SCNS03 & PR-SCNS04 disposable shoe covers are made of a polypropylene fabric and are 100% latex and silicone free. The high-quality fabric and strong welds set them apart from other shoe covers. They are constructed with tripleelastic runs that allow a smooth wrap around the foot in order to prevent tripping. They are excellent for transitioning between "clean" and "unclean" areas.

Features

- Elastic band for secure
- Triple-elastic runs for s
- Lightweight but durable
- Slip-resistant pattern
- Universal and 2XL Sizes
- Easy Open Bag
- Fluid resistant
- Color: Blue



PR-SCSB03 & PR-SCSB02

The PR-SCSB03 & PR-SCSB02 disposable shoe covers are made of a proprietary fabric that has higher wear resistance, fluid resistance and antiskid properties for better foot grip. They are tested for linting compliance for cleanrooms up to class ISO 5 (FED STD 209E Class 100).

Features

- Elasticband for secure fit
- Waterproof
- Slip-resistant
- Universal and XL Sizes
- Easy Open Bag
- Very low Linting
- Color: Light Blue



ICP-BCNC0003

The ICP-BCNC0003 is an impervious disposable boot cover made of a polypropylene fabric with two conductive strips on the bottom. It is 100% latex.

Features

- Elastic band at knee
- Tie at ankle
- Waterproof
- Non-Conductive Strips
- Universal
- Color: Blue



Protexer | Disposable Personal Protection Equipment

Bouffant Caps, Sleeve Covers and Beard Covers

Description

Our bouffant caps are available in light weight single ply polypropylene or a heavy weight spunbound polypropylene. Secured with a latex free elastic band and available in multiple sizes.

Bouffant Caps:

ICP-CPPSB0021 Fluid Resistant Bouffant Cap, 21", 1000/Case

ICP-CPPSB0024 Fluid Resistant Bouffant Cap, 24", 1000/Case

ICP-CPPSB0028 Fluid Resistant Bouffant Cap, 28", 500/Case

ICP-CPXP24 Extra Protection Bouffant Cap, 24", 500/Case



Description

Sleeve covers are constructed of a PP/PE material creating a high fluid resistant barrier.

Sleeve Cover:

PR-SCW-500 Sleeve Cover, PP/PE, White, 500/Case



Description

Beard Covers made of a spunbound polypropylene, stays in place with an elastic band

Beard Cover:

PR-BCB-500 Beard Cover, Blue, 500/Case



Protexer | Stretcher Sheets

Ultra-protective Stretcher Sheets

These disposable fluid-resistant, fitted sheets provide protection from patient secretions and body fluids. The soft-touch fabric is ultra comfortable for patients. Convenient, easy to change and sized to fit most gurneys, these fitted sheets can save valuable nursing time.


- Ultra-protective
- Minimizes risk of cross-contamination
- Fluid-repellant SPP fabric
- Elastic ends with deep corner pockets
- Tear-resistant
- Great backup for laundered sheets
- Reduces laundry/linen costs
- Latex-free

ICP-9001GUR

96" x 37" white SPP, 1 pc/bag;
100 pcs/case

ICP-9003GUR

108" x 51" white SPP, 1 pc/bag;
100 pcs/case



TEAM Medical, a division of TEAM Technologies, brings together a portfolio of solutions for the healthcare industry with a variety of products ranging from PPE to specialized devices. Our TEAM of experts offers end-to-end solutions for manufacturing, supply chain and distribution needs as well as products used in healthcare settings ranging from the clinic to commercial applications.

Our Brands/Products:

- **Rapid Refresh Privacy Curtains**
- **Protexer PPE Apparel and disposable products**
 - Gowns, Scrubs, Shoe Covers, Bouffant Caps, Lab Coats, Sleeve Covers, Beard Covers, Stretcher Sheets
- **Protexer PPE Organizers**
 - Customizable Configurations, Glove Boxes, Enclosed Cabinets, PPE Dispensers
- **Protexer Sticky Mats**
- **BootieButler Shoe Covers**
- **DoseLogix**
 - Topi-CLICK
- **Pure Disinfectant**
- **ICP Medical Keyboards**

For more information contact:

TEAM Technologies

5949 Commerce Blvd
Morristown, TN 37814
toll-free: **1-423-587-2199**
sales@teamtech.com
www.teamtech.com

ICP Medical

13720 Rider Trail North
St. Louis, MO 63045
toll-free: **1-800-405-3044**
info.ICPMedical@teamtech.com
www.icpmedical.com

BootieButler

1616 Park 370 Ct.
St. Louis MO 63042
toll-free: **1-800-710-9863**
info.BootieButler@teamtech.com
www.bootiebutler.com

DoseLogix

300 Parkbrook Place, Suite 140
Woodstock, GA 30189
toll-free: **1-877-870-8448**
info.DoseLogix@teamtech.com
www.doselogix.com

www.teamtech.com

